Osmometer Standards

REF	
30.9.0100	Zehn 1-mL Ampullen à 100 mOsmol/kg NaCl/H ₂ O
30.9.0290	Zehn 1-mL Ampullen à 290 mOsmol/kg NaCl/H ₂ O
30.9.0020	Zehn 1-mL Ampullen à 300 mOsmol/kg NaCl/H ₂ O
30.9.0500	Zehn 1-mL Ampullen à 500 mOsmol/kg NaCl/H ₂ O
30.9.0850	Zehn 1-mL Ampullen à 850 mOsmol/kg NaCl/H ₂ O
30.9.2000	Zehn 1-mL Ampullen à 2000 mOsmol/kg NaCl/H ₂ O

IVD

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IVD

INTENDED USE

Gonotec calibration standards are used to calibrate and/or check the performance of the osmometer by qualified personnel.

This product meets the requirements of EEC Directive 98/79 relating to in-vitro diagnostics

INGREDIENTS

The standards are prepared from NaCl (clinical standards) and water for injection purpose. The salt meets the requirements of EUR, PH., the US PH. and the JP XIV. The ampoules are sterilized in auto claves. Concentration data can be found on the ampoule and the outer package.

WARNING

"One Point Cut" (OPC) – ampoules made of colourless glass are used with a blue marked breakpoint. At that point the ampoule can manually be opened by snapping off. All safety measures regarding the use of glass have to be observed (building of splinter, break of glass, ...).

CLASS OF DANGER

Natriumchloride The calibration- and reference material is not a CAS-NR.:7647-14-5 (NaCl) dangerous substance according to EC Directive 67/548/EEC.

Wasser

CAS-NR.:7732-18-5 (H₂O)

Poison-class (CH)

F (free)

Water-endangering class

ne (no endangering)

Storage-class

LGK 10-13: (other liquids and solids)

R-phrase

None

S-phrase

S 1/2 To be stored locked and kept out of reach of children.

S 24/25 Avoid contact with eyes and skin.

S 26 After contact rinse eyes thoroughly with water, see an eye specialist.

STORAGE AND STABILITY

	Storage	Stability
Unopened	5 - 45°C (41 - 113°F)	Refer to expiry date on the ampoule or outer package
Opened	22°C (71,6°F)	max ½ hour

- Calibration readings are controlled by repeated calibrations. Increasing osmolality is a sign for exceeding stability.
- Do not freeze opened ampoules.
- Do not mix remaining solvents of opened ampoules. No reuse.

LIMITATIONS

Reuse or use of expired calibrators, adverse shipping and/or storage conditions can lead to erroneous results as well as the use of defective pipettes, dirty syringes and/or sample vessels or erroneous use of the samples.

CALIBRATION

The sample volume and the calibration process can be learned from the corresponding User's Guide of the device.

PREPARATION

- Carefully snap off the top of the ampoule.
- Obtain samples with calibrated pipette with clean tip.
- Refer to User's Guide for more information about sample size and testing procedure.

QUALITY CONTROL

Gonotec calibration and reference solutions are required for monitoring sample measurements as well as the quality of the calibration curve. Upon evaluation of control values and the determination of admissible ranges appropriate statistical methods should be used. After calibration reference measurements should be made, e.g. with other solutions with known osmolality to verify the acceptance of the calibration. At least one reference material of known osmolality should be used as control. Deviations to the expected result may be a sign for erroneous use of the device, the calibration and reference material, for contaminated syringes or vessels or a malfunction of the device.

CHARACTERISTICS OF THE CALIBRATION

The table below lists the nominal value and the corresponding expected ranges for the osmometer standards, assuming nominal instrument performance.

Expected Value	Expected range (bias)
100 mOsmol/kg NaCl/H ₂ O	98 – 102 mOsmol/kg
290 mOsmol/kg NaCl/H ₂ O	288 – 292 mOsmol/kg
300 mOsmol/kg NaCl/H ₂ O	298 – 302 mOsmol/kg
500 mOsmol/kg NaCl/H ₂ O	497.5 – 502.5 mOsmol/kg
850 mOsmol/kg NaCl/H ₂ O	845.7 – 854.3 mOsmol/kg
2000 mOsmol/kg NaCl/H ₂ O	1990 – 2010 mOsmol/kg

DISPOSING OF MATERIALS

Handle this product according to established good laboratory practise (GLP). We recommend transferring the solution in containers for neutral aqueous solutions. Observe local requirements for disposal of chemicals and rest material. Never dispose remainders in house rubbish/sink/WC.

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CAL

Osmometer Standard

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Achtung, Dokumentation beachten!
Attention, refer to documentation!

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FOR SALES AND SERVICE:
Please Contact Your Gonotec
Distributor

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Manufactured and Printed in Germany
As at December 2010

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Gonotec GmbH