

# 吸头适配表

Pipette Brand		Axygen Axypet (Single channel)						Axygen Axypet (Multichannel)				Nichiyu/Oxford Benchmate			Biohit Proline (Single channel)				
Model / Range		0.1-2µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	1-10µl	5-50µl	20-200µl	50-300µl	0.5-10µl	2-20µl	10-50µl	0.2-10µl	5-100µl	10-250µl	10-500µl	50-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓					✓				✓			✓				
	21-0010	✓	✓					✓				✓			✓				
1-200ul	20-0200			✓	✓	✓	✓		✓	✓	✓		✓	✓		✓	✓		
	21-0200			✓	✓	✓	✓		✓	✓	✓		✓	✓		✓	✓		
100-1000ul	20-1000																	✓	✓
	21-1000																	✓	✓
100-1250ul	20-1250																	✓	✓
	20-1251																	✓	✓
	21-1250																	✓	✓

Pipette Brand		Biohit Proline			Biohit e-line (Single channel)				Biohit e-line (Multi)			Biohit m-line (Single channel)			Biohit m-line				
Model / Range		0.2-10µl	5-100µl	25-250µl	1-10µl	5-120µl	10-300µl	50-1000µl	0.2-10µl	5-120µl	10-300µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	0.5-10µl	5-100µl	30-300µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓			✓				✓									✓	
	21-0010	✓			✓				✓									✓	
1-200ul	20-0200		✓	✓		✓	✓			✓			✓	✓				✓	✓
	21-0200		✓	✓		✓	✓			✓			✓	✓				✓	✓
100-1000ul	20-1000							✓								✓			
	21-1000							✓								✓			
100-1250ul	20-1250							✓								✓			
	20-1251							✓								✓			
	21-1250							✓								✓			

Pipette Brand		Capp Aero (Multichannel)							BrandTech Transferpette (Single channel)										
Model / Range		0.2-2µl	0.5-10µl	5-50µl	10-100µl	25-200µl	50-300µl	100-1000µl	0.1-1µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	25-250µl	100-1000µl	20-200µl	5-100µl	30-300µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓						✓	✓									
	21-0010	✓	✓						✓	✓									
1-200ul	20-0200			✓	✓	✓	✓				✓	✓	✓	✓			✓		
	21-0200			✓	✓	✓	✓				✓	✓	✓	✓			✓		
100-1000ul	20-1000							✓							✓				
	21-1000							✓							✓				
100-1250ul	20-1250							✓							✓				
	20-1251							✓							✓				
	21-1250							✓							✓				✓

# 吸头适配表

Pipette Brand		BrandTech Transferpette (Multichaanel)						Rainin EDP3 (Single channel)						Rainin Latch Mode					
Model / Range		0.5-10µl	2-20µl	2.5-25µl	5-50µl	10-100µl	20-200µl	30-300µl	E-MIC-10, 0.5-10µl	E2-25, 2.5-25µl	E2-100, 10-100µl	E2-250, 25-250µl	E2-1000, 100-1000µl	R2, 0.1-2µl	R10, 0.2-10µl	R20, 2-20µl	R100, 10-100µl	R200, 20-200µl	R1000, 100-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓						✓					✓	✓				
	21-0010	✓	✓						✓					✓	✓				
1-200ul	20-0200			✓	✓	✓	✓		✓	✓	✓				✓	✓	✓	✓	
	21-0200			✓	✓	✓	✓		✓	✓	✓				✓	✓	✓	✓	
100-1000ul	20-1000							✓					✓						✓
	21-1000							✓					✓						✓
100-1250ul	20-1250							✓					✓						✓
	20-1251							✓					✓						✓
	21-1250							✓					✓						✓

Pipette Brand		Labnet Labpette (Single chaanel)					Labnet Labpette (Multichaanel)				Labnet Biopette E (Single chaanel)			Labnet Biopette E					
Model / Range		0.1-2µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	0.5-10µl	2-20µl	10-200µl	100-1000µl	0.5-10µl	2-20µl	10-200µl	100-1000µl
Volume	Product No.	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓					✓				✓	✓		✓	✓			
	21-0010	✓	✓					✓				✓	✓		✓	✓			
1-200ul	20-0200			✓	✓	✓		✓	✓	✓			✓					✓	
	21-0200			✓	✓	✓		✓	✓	✓			✓					✓	
100-1000ul	20-1000						✓							✓					✓
	21-1000						✓							✓					✓
100-1250ul	20-1250						✓							✓					✓
	20-1251						✓							✓					✓
	21-1250						✓							✓					✓

Pipette Brand		Eppendorf Electronic				Socorex Calibra/Accura (Single chaanel)						VWR AND VWR Signature (Multichannel)				Rainin EDP3 (Multichannel)				
Model / Range		0.5-10µl	5-100µl	2-200µl	50-1000µl	0.1-2µl	0.2-2µl	0.5-10µl	1-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	1-10µl	5-50µl	20-200µl	50-300µl	EP-M8-10, 0.5-10µl	EP-M8250, 10-250µl	P200-M8, 20-200µl
Volume	Product No.	Universal Fit Filter Tips																		
0.5-10ul	20-0010	✓				✓	✓	✓	✓					✓				✓		
	21-0010	✓				✓	✓	✓	✓					✓				✓		
1-200ul	20-0200		✓	✓					✓	✓	✓			✓	✓	✓			✓	✓
	21-0200		✓	✓					✓	✓	✓			✓	✓	✓			✓	✓
100-1000ul	20-1000				✓								✓							✓
	21-1000				✓								✓							✓
100-1250ul	20-1250				✓								✓							✓
	20-1251				✓								✓							✓
	21-1250				✓								✓							✓

# 吸头适配表

Pipette Brand		VWR AND VWR Signature (Single channel)							Socorex Calibra/Accura (Multichannel)					EppendorfReference (single channel)					
Model / Range		0.1-2µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	50-250µl	100-1000µl	1-10µl	5-50µl	10-100µl	20-200µl	50-200µl	0.1-10µl	2-20µl	10-100µl	20-200µl	100-1000µl
Volume	Product No	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓							✓					✓				
	21-0010	✓	✓							✓					✓				
1-200ul	20-0200			✓	✓	✓	✓	✓		✓	✓	✓	✓		✓	✓	✓	✓	
	21-0200			✓	✓	✓	✓	✓		✓	✓	✓	✓		✓	✓	✓	✓	
100-1000ul	20-1000								✓						✓	✓	✓	✓	
	21-1000								✓						✓	✓	✓	✓	
100-1250ul	20-1250								✓										✓
	20-1251								✓										✓
	21-1250								✓										✓

Pipette Brand		Eppendorf Research					Eppendorf Research (multichannel)					Thermo Scientif Finnpiquette (single channel)							
Model / Range		0.1-2.5µl	0.5-10µl	2-20µl	10-100µl	25-200µl	100-1000µl	0.5-10µl	5-100µl	50-100µl	20-300µl	30-300µl	50-1200µl	0.2-10µl	0.5-10µl	5-40µl	40-200µl	100-1000µl	200-1000µl
Volume	Product No	Universal Fit Filter Tips																	
0.5-10ul	20-0010	✓	✓					✓						✓	✓				
	21-0010	✓	✓					✓						✓	✓				
1-200ul	20-0200				✓	✓			✓	✓	✓	✓			✓	✓	✓	✓	
	21-0200				✓	✓			✓	✓	✓	✓			✓	✓	✓	✓	
100-1000ul	20-1000				✓	✓			✓	✓	✓	✓			✓	✓	✓	✓	
	21-1000						✓					✓						✓	✓
100-1250ul	20-1250						✓					✓						✓	✓
	20-1251						✓					✓						✓	✓
	21-1250						✓					✓						✓	✓

Pipette Brand		Thermo Scientif Finnpiquette		Gilson Pipetma (single channel)					
Model / Range		5-50µl	50-300µl	P20, 1-2µl	P10, 0.2-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl
Volume	Product No	Universal Fit Filter Tips							
0.5-10ul	20-0010			✓	✓				
	21-0010			✓	✓				
1-200ul	20-0200	✓	✓			✓	✓	✓	
	21-0200	✓	✓			✓	✓	✓	
100-1000ul	20-1000	✓	✓			✓	✓	✓	
	21-1000								✓
100-1250ul	20-1250								✓
	20-1251								✓
	21-1250								✓

Pipette Brand		Eppendorf Research						Gilson Pipetman (PG)						Gilson Pipetman (PN)						Gilson Pipetman (PL)					
Model / Range		0.1-2.5µl	0.5-10µl	2-20µl	10-100µl	20-200µl	100-1000µl	P2, 0.2-2µl	P10, 1-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl	P2, 0.2-2µl	P10, 1-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl	P2, 0.2-2µl	P10, 1-10µl	P20, 2-20µl	P100, 10-100µl	P200, 20-200µl	P1000, 100-1000µl
Volume	Product No	Universal Fit Filter Tips																							
0.5-10ul	22-0010S	✓	✓					✓	✓					✓	✓					✓	✓				
	22-0011S	✓	✓					✓	✓					✓	✓					✓	✓				
1-20ul	22-0020S			✓						✓													✓		
	22-0021S			✓						✓													✓		
1-50ul	22-0050S				✓						✓	✓					✓	✓					✓	✓	
	22-0051S				✓						✓	✓					✓	✓					✓	✓	
1-100ul	22-0100S				✓						✓	✓					✓	✓					✓	✓	
	22-0101S				✓						✓	✓					✓	✓					✓	✓	
1-300ul	22-0200S					✓						✓						✓	✓					✓	
	22-0201S					✓						✓						✓	✓					✓	
100-1000ul	22-1000S						✓					✓						✓	✓					✓	
	22-1001S						✓					✓						✓	✓					✓	

Pipette Brand		Thermo Finnpipette (F1)										Thermo Finnpipette										Thermo Finnpipette													
Model / Range		0.2-2µl	1-10µl	1-10µl(µlTip)	2-20µl	2-20µl(µlTip)	5-50µl	5-50µl(µlTip)	10-100µl	20-200µl	30-300µl	100-1000µl	0.2-2µl(µlTip)	1-10µl	1-10µl(µlTip)	2-20µl	2-20µl(µlTip)	5-50µl	5-50µl(µlTip)	10-100µl	20-200µl	100-1000µl	0.2-2µl	1-10µl	1-10µl(µlTip)	2-20µl	2-20µl(µlTip)	5-50µl	5-50µl(µlTip)	10-100µl	20-200µl	100-1000µl			
Volume	Product No	Universal Fit Filter Tips																																	
0.5-10ul	22-0010S	✓	✓											✓									✓	✓											
	22-0011S	✓	✓											✓									✓	✓											
1-20ul	22-0020S			✓												✓										✓									
	22-0021S			✓												✓										✓									
1-50ul	22-0050S				✓			✓	✓							✓	✓									✓			✓	✓					
	22-0051S				✓			✓	✓							✓	✓									✓			✓	✓					
1-100ul	22-0100S				✓			✓	✓							✓	✓									✓			✓	✓					
	22-0101S				✓			✓	✓							✓	✓									✓			✓	✓					
1-300ul	22-0200S					✓					✓																		✓	✓					
	22-0201S					✓					✓																		✓	✓					
100-1000ul	22-1000S						✓					✓												✓						✓					
	22-1001S						✓					✓												✓						✓					

# 吸头适配表

Pipette Brand		Thermo Finnpipette (F1)								Thermo Finnpipette								Thermo Finnpipette									
Model / Range		0.1-2.5µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	50-300µl	100-1000µl	0.1-3µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	30-300µl	100-1000µl	0.1-1µl	0.1-2.5µl	0.5-10µl	2-20µl	5-50µl	10-100µl	20-200µl	30-300µl	100-1000µl	
Volume	Product No.	Universal Fit Filter Tips																									
0.5-10ul	22-0010S	✓	✓							✓	✓							✓	✓	✓							
	22-0011S	✓	✓							✓	✓							✓	✓	✓							
1-20ul	22-0020S			✓							✓									✓							
	22-0021S			✓							✓									✓							
1-50ul	22-0050S				✓	✓	✓	✓				✓	✓	✓	✓						✓	✓	✓	✓			
	22-0051S				✓	✓	✓	✓				✓	✓	✓	✓						✓	✓	✓	✓			
1-100ul	22-0100S					✓	✓	✓					✓	✓	✓							✓	✓	✓	✓		
	22-0101S					✓	✓	✓					✓	✓	✓								✓	✓	✓	✓	
1-300ul	22-0200S						✓	✓							✓	✓								✓	✓		
	22-0201S						✓	✓							✓	✓									✓	✓	
100-1000ul	22-1000S								✓								✓										✓
	22-1001S								✓								✓										✓

Pipette Brand		Eppendorf Research Plus							Gilson Pipetman (PG)					Gilson Pipetman (PN)					Gilson Pipetman (PL)								
Model / Range		0.1-2.5ul	0.5-10ul	2-20ul	10-100ul	20-200ul	100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul	P2, 0.2-2ul	P10, 1-10ul	P20, 2-20ul	P100, 10-100ul	P200, 20-200ul	P1000, 100-1000ul		
Volume	Product No.	Universal Fit Filter Tips																									
0.5-10ul	23-0010S	✓	✓					✓	✓					✓	✓					✓	✓						
	23-0011S	✓	✓					✓	✓					✓	✓					✓	✓						
1-20ul	23-0020S			✓					✓						✓							✓					
	23-0021S			✓					✓						✓							✓					
1-50ul	23-0050S				✓					✓	✓					✓	✓					✓	✓	✓	✓		
	23-0051S				✓					✓	✓					✓	✓					✓	✓	✓	✓		
1-100ul	23-0100S				✓					✓	✓					✓	✓					✓	✓	✓	✓		
	23-0101S				✓					✓	✓					✓	✓					✓	✓	✓	✓		
1-300ul	23-0200S					✓						✓						✓						✓	✓		
	23-0201S					✓						✓						✓							✓	✓	
100-1000ul	23-1000S						✓						✓						✓								✓
	23-1001S						✓						✓						✓								✓

